





High-Definition Wireless Ultrasound for Regional Anesthesia

Clear, Rapid Guidance

With high-definition imaging comparable to traditional ultrasound systems, Clarius is the ultimate wireless ultrasound scanner for safe, longer-lasting regional blocks and post-operative follow-up anywhere.



Visualize Your Needle in High-Definition for Safe Procedures

Optimize patient comfort and minimize complications. Wireless and pocket-sized, Clarius delivers the high-definition imaging and performance of traditional systems in a highly affordable ultra-mobile scanner. With 8 times the processing power of most handhelds, you get clear visualization of your needle and target anatomy for safe, accurate procedures, every time.

- Use ultrasound to guide a supraclavicular block in the pre-op area
- Visualize an intra-articular injection in the operating room
- Confirm correct positioning of a continuous catheter
- Run to the post-anesthetic care unit for a pectoralis nerve block
- Visualize your needle, nerve, and vessels for every procedure



**Clarius L7 HD3
Linear Scanner**

**Frequency: 4-13 MHz
Max Depth: 11 cm**




**Clarius L15 HD3
High-Freq. Scanner**

**Frequency: 5-15 MHz
Max Depth: 7 cm**




**Clarius C3 HD3
Curvilinear Scanner**


**Frequency: 2-6 MHz
Max Depth: 40 cm**



Axillary Nerve




Supraclavicular Brachial Plexus Block




Supraclav Brachial Plexus



Infraclavicular Plexus



Axillary Nerve Injection



Adductor Canal Block



Wireless ultrasound enables me to focus more on my patient without having wires get in the way of my nerve block. The third generation Clarius is amazing. The image quality is superior. It connects seamlessly. I mean, it's light, it's small, and it's easier to handle when performing a procedure. It's cheaper, more portable and has greater flexibility than the traditional ultrasound machines."

Dr. David Rosenblum, MD - Anesthesiologist

Image Clarity You Can Trust for Regional Anesthesia



High-Definition

8 beamformers, 192 piezoelectric elements, and AI deliver high performance and resolution.



Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.



Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



Ultra-Portable

Clarius fits in your pocket to guide safe regional blocks anywhere.



Specialized

AI-powered with specialized MSK presets for anesthesia.




Exam Management

Seamlessly transfer saved images and videos to PACS servers through the current DICOM web standard.




The beauty of the wireless ultrasound is it's so mobile, lightweight, and it's so easy to use. It improves your speed of being able to do the block, and it actually improves the duration of the block because you're getting more local around the target nerve versus doing it blind."

Dr. Gregory Hickman - Anesthesiologist



Advanced MSK Package

The Advanced MSK Package offers more flexibility for users who need additional customizations for examinations and procedures for various joint spaces and MSK related anatomy.



Dedicated Presets

Optimize imaging with presets for shoulder, elbow, hand/wrist, knee, foot/ankle, hip, spine, and more.



Needle Enhance

Enhance the visualization of needles when performing guided procedures.



Specialized Measurements

Save time and improve reporting by having pre-labeled measurements at your disposal.

Additional Features

- DICOM
- Elastography
- Split Screen Mode

Clarius Membership

The Advanced MSK Package is included with your Clarius Membership.

Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Education
- ✓ Advanced MSK Package

One-Time Purchase

You may also purchase your scanner outright with the Advanced MSK Package à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- 💰 Advanced MSK Package

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

Clarius Mobile Health

2985 Virtual Way #130,
Vancouver, BC, V5M 4X7
Canada
hello@clarius.com
Phone 1-778-800-9975