



High-Definition Wireless Ultrasound for Veterinarians

Fast Confident Diagnosis

Find answers under the fur with high-definition wireless ultrasound. Designed for veterinarians, Clarius HD3 Vet wireless scanners bring AI-powered imaging with dedicated animal presets to your iOS or Android device for a confident diagnosis anywhere.



Speed Diagnosis to Expedite Treatment for Your Furry Patients

With high-definition imaging and automated settings, Clarius HD3 Vet wireless scanners are the ultimate POCUS solution for animals small, medium, and large. Affordable, AI-powered, and with dedicated animal presets, Clarius is ultra-portable and streamlines workflows for a confident diagnosis anywhere. Make every room an exam room, scanning felines by the cage side and large canines on the floor, to get all your furry patients back in action faster!

- Perform an abdominal FAST exam to identify free fluid
- Take accurate measurements of masses, sludge, and free fluid
- Scan the adrenal glands for hyperplasia, nodules, and masses
- Scan the lung and pleural space for pathology of cardiac origin
- Assess the heart to determine if furosemide is indicated
- Perform an equine exam of the superficial digital flexor tendon



Clarius C7 HD3 VET
Small - Medium

Frequency: 3-10 MHz
Max Depth: 15 cm



Clarius C3 HD3 VET
Medium - Large

Frequency: 2-6 MHz
Max Depth: 40 cm




Clarius L7 HD3 VET
Animal MSK


Frequency: 4-13 MHz
Max Depth: 7 cm




Canine Cardiac (C7)




Dog Kidney (C7)




Equine Tendon (L7)



Feline Bladder (C7)



Canine Pleural Effusion (C7)



Equine Facet Joint (C3)



I have certainly recommended it to colleagues who are considering a handheld ultrasound machine. The Clarius is the best one that I've tried. The image quality is fantastic for the price. If you are new to ultrasound, it is very easy to get started with because the AI sorts out a lot of the knobs and functions that you would have to learn and adjust manually on another ultrasound machine."

Dr. Camilla Edwards, DVM, CertAVP, MRCVS - Ultrasound Educator

Specialized Scanners for High-Resolution Ultrasound



Clarius Cloud

Manage your exams and build reports anywhere, anytime with unlimited Clarius Cloud storage.



Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



High-Definition

Latest Octal Beamforming and AI technology provides best-in-class image quality.



Specialized

AI-powered with specialized workflows for veterinarians.




Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!



Wireless

Steadfast connection. Total freedom. Without cords to startle animals, cleaning is also simple and fast.




Clarius ultrasound has changed the way we look at animals and enabled us to provide better care. Animals don't like being handled and having a lot of equipment around. So, having a small system with an iPad that they're already familiar with helps us to get an exam done quickly. Since it's so compact, we can carry it with us anywhere. It used to be difficult to carry our trolley system up and down the stairs."

Dr. Luke Ross, DMV - Veterinarian




Advanced Veterinary Package

The Advanced Veterinary Package offers more flexibility for veterinarians who want advanced workflows for various animal examinations. Artificial Intelligence fine tunes imaging to visualize finer details and provide better care with specialized presets.



Advanced Presets

Take veterinary scanning to the next level by gaining access to finely tuned, categorized presets.



Specialized Measurements

Improve workflows and generate detailed reports with animal-driven measurements and calculations.

Additional Features

- DICOM
- Elastography
- Needle Enhance
- Split Screen Mode
- Pulsed Wave Doppler

Clarius Membership

The Advanced Veterinary Package is included with your Clarius Membership.

Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Education
- ✓ Advanced Veterinary Package

One-Time Purchase

You may also purchase your scanner outright with the Advanced Veterinary Package à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✶ Advanced Veterinary Package

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

Clarius Mobile Health

2985 Virtual Way #130,
Vancouver, BC, V5M 4X7
Canada
hello@clarius.com
Phone 1-778-800-9975